

35MM ENDOCUTTERS AND RELOADS







EES DEVICE	CODE	DESCRIPTION	ARTICUL.	RELOAD	MAX. FIRINGS	NON-EES DEVICE
ETS 35MM	TSB35 ●	Linear Cutter w/ 6 Row Blue Reload	NO	TR35B ●	8	ENDO GIA™ UNIVERSAL
	TSW35 ○	Linear Cutter w/ 6 Row White Reload		TR35W ○		
FLEX 35MM	ATB35 ●	Linear Cutter w/ 6 Row Blue Reload	YES			
	ATW35 ○	Linear Cutter w/ 6 Row White Reload				

ETS 35mm and FLEX 35mm reloads fit down 12mm trocar.

35MM ENDOCUTTER RELOADS



EES DEVICE	CODE	DESCRIPTION	OPEN STAPLE HT.	CLOSED STAPLE HT.	NON-EES DEVICE
ETS 35MM	TR35W ○	White, 6 Row, Vascular/Thin Tissue	 2.5mm	 1.0mm	ENDO GIA™ 30mm Reloads
	TR35B ●	Blue, 6 Row, Regular Tissue	 3.5mm	 1.5mm	

LAP
OPEN



45MM ENDOCUTTERS AND RELOADS













EES DEVICE	CODE	DESCRIPTION	ARTICUL.	RELOAD	MAX. FIRINGS	NON-EES DEVICE
ETS 45MM	ETS45	Regular Linear Cutter	NO	6R45M ●	8	ENDO GIA™ UNIVERSAL
				TR45W ○		
FLEX 45MM	ATS45	FLEX Regular Linear Cutter	YES	TR45B ●		
				6R45B ●		
				TR45G ●		
				CTS45	FLEX Compact Linear Cutter	
				ATS45NK	FLEX Regular Linear Cutter, No Knife	
CTS45NK	FLEX Compact Linear Cutter, No Knife					
	LONG45A	FLEX Long Linear Cutter				

ETS endocutters do not contain a reload. ETS 45mm fits down 12mm trocar.

45MM ENDOCUTTER RELOADS



EES DEVICE	CODE	DESCRIPTION	OPEN STAPLE HT.	CLOSED STAPLE HT.	NON-EES DEVICE	
ETS 45MM	6R45M ●	Gray, 6 Row, Mesentery/Thin Tissue	 2.0mm	 .85mm	ENDO GIA™ 45mm Reloads ^{1, 2, 3}	
	TR45W ○	White, 6 Row, Vascular/Thin Tissue	 2.5mm	 1.0mm		Roticulating or Straight
	TR45B ●	Blue, 4 Row, Regular Tissue	 3.5mm	 1.5mm		
	6R45B ●	Blue, 6 Row, Regular Tissue	 3.5mm	 1.5mm		
	TR45G ●	Green, 4 Row, Thick Tissue	 4.1mm	 2.0mm		

ETS 45mm reloads fit down 12mm trocar.

- 1. ENDO GIA™ gray closed staple height .75mm.
- 2. ENDO GIA™ green reload requires a 15mm trocar.
- 3. Green open staple height is 4.8mm. Same closed staple height.

LAP
OPEN

60MM ENDOCUTTERS AND RELOADS



EES DEVICE	CODE	DESCRIPTION	ARTICUL.	RELOAD	MAX. FIRINGS	NON-EES DEVICE
ETS 60MM	LTS60A	FLEX ENDO Linear Cutter	YES	TR60W ○ TR60B ● TR60G ●	8	ENDO GIA™ UNIVERSAL

ETS endocutters do not contain a reload. ETS 60mm and ETS 60mm reloads fit down 12mm trocar.

60MM ENDOCUTTER RELOADS



EES DEVICE	CODE	DESCRIPTION	OPEN STAPLE HT.	CLOSED STAPLE HT.	NON-EES DEVICE
ETS 60MM	TR60W ○	White, 6 Row, Vascular/Thin Tissue	 2.5mm	 1.0mm	ENDO GIA™ 60mm Reloads ^{1,2}
	TR60B ●	Blue, 6 Row, Regular Tissue	 3.5mm	 1.5mm	Roticulating or Straight
	TR60G ●	Green, 4 Row, Thick Tissue	 4.1mm	 2.0mm	

1. ENDO GIA™ green reload requires a 15mm trocar. 2. Green open staple height is 4.8mm. Same closed staple height.