

LINEAR STAPLERS AND RELOADS

PROXIMATE®



EES DEVICE	CODE	DESCRIPTION	RELOAD	MAX. FIRINGS	NON-EES DEVICE
PROXIMATE® 30MM	TX30V	RL Vascular, Secures Tissue Thickness of 1.0mm	XR30V	8	DST Series™ TA™ 30mm
	TX30B	RL, Linear Stapler Secures Tissue Thickness of 1.5mm	XR30B XR30G	8	
	TX30G	RL, Linear Stapler Secures Tissue Thickness of 2.0mm			
COMPRESSION DIAL	TLV30	RL, Titanium Vascular, Secures Tissue Thickness of 1.0mm	TRV30	4	
	TL30	RL Titanium, Secures Tissue Thickness of 1.0 - 2.5mm	TR30	4	
	TLH30	RL Titanium, Heavy Wire Secures Tissue Thickness of 1.0 - 2.5mm	TRH30	4	
60MM	TX60B	RL, Linear Stapler Secures Tissue Thickness of 1.5mm	XR60B XR60G	8	DST Series™ TA™ 60mm
	TX60G	RL, Linear Stapler Secures Tissue Thickness of 2.0mm			
COMPRESSION DIAL	TL60	RL Titanium, Secures Tissue Thickness of 1.0 - 2.5mm	TR60	4	
	TLH60	RL Titanium, Heavy Wire Secures Tissue Thickness of 1.5 - 2.5mm	TRH60	4	
90MM	TL90	RL Titanium, Secures Tissue Thickness of 1.5 - 2.5mm	TR90	4	DST Series™ TA™ 90mm
	TLH90	RL Titanium, Heavy Wire Secures Tissue Thickness of 1.0 - 2.5mm	TRH90	4	

LAP

OPEN